

# Quantikine<sup>®</sup>

## Human BMP-4 Immunoassay

Catalog Number DBP400

**For the quantitative determination of human bone morphogenetic protein 4 (BMP-4) concentrations in bone tissue extracts and cell culture supernates.**

*This package insert must be read in its entirety before using this product.*

**FOR RESEARCH USE ONLY.  
NOT FOR USE IN DIAGNOSTIC PROCEDURES.**

# TABLE OF CONTENTS

Contents	Page
INTRODUCTION	2
PRINCIPLE OF THE ASSAY . . . . .	3
LIMITATIONS OF THE PROCEDURE	3
REAGENTS . . . . .	3
STORAGE	4
OTHER SUPPLIES REQUIRED . . . . .	4
PRECAUTION	4
SAMPLE COLLECTION AND STORAGE . . . . .	5
REAGENT PREPARATION	5
ASSAY PROCEDURE . . . . .	6
ASSAY PROCEDURE SUMMARY	7
CALCULATION OF RESULTS . . . . .	8
TYPICAL DATA	8
TECHNICAL HINTS . . . . .	9
PRECISION	9
RECOVERY . . . . .	9
SENSITIVITY	10
CALIBRATION . . . . .	10
LINEARITY	10
SPECIFICITY . . . . .	10
REFERENCES	11

## MANUFACTURED AND DISTRIBUTED BY:

R&D Systems, Inc.	TELEPHONE:	(800) 343-7475
614 McKinley Place N.E.		(612) 379-2956
Minneapolis, MN 55413	FAX:	(612) 379-6580
United States of America	E-MAIL:	info@RnDSystems.com

## DISTRIBUTED BY:

R&D Systems Europe	TELEPHONE:	(0)1235 529449
19 Barton Lane	FAX:	(0)1235 533420
Abingdon Science Park	E-MAIL:	info@RnDSystems.co.uk
Abingdon, Oxon OX14 3NB		
United Kingdom		

R&D Systems GmbH	FREEPHONE:	(0)800 909 4455
Borsigstrasse 7	TELEPHONE:	(0)6122 90980
65205 Wiesbaden-Nordenstadt	FAX:	(0)6122 909819
Germany	E-MAIL:	infogmbh@RnDSystems.co.uk

## INTRODUCTION

Bone morphogenetic protein-4 (BMP-4, previously known as BMP-2b) is a member of the transforming growth factor beta (TGF- $\beta$ ) superfamily. BMPs were originally identified as protein regulators of cartilage and bone formation. They have also been implicated in embryogenesis and morphogenesis of various tissues and organs. They can regulate growth, differentiation, chemotaxis and apoptosis of a variety of cell types, including mesenchymal, epithelial, hematopoietic and neuronal cells. BMP-4 plays an important role in the onset of endochondral bone formation, dorsal/ventral patterning and has also been implicated in the commitment of embryonic mesodermal cells to a hematopoietic fate in a number of systems (for reviews, see references 1 - 5).

Each BMP is synthesized as a precursor peptide, processed to a mature form, and subsequently secreted as a dimer. Although homodimers are considered the standard form, there are natural heterodimers with equal, if not increased, bioactivity (6, 7). BMP-4 is a 408 amino acid (aa) prepropeptide composed of a 19 aa signal sequence, a 273 aa pro-region, and a 116 aa mature segment (8). Both the pro-region and mature segment contain two potential N-linked glycosylation sites. The mature region contains seven highly conserved cysteine aa residues that form the characteristic cysteine knot found within TGF- $\beta$  superfamily members. The mature regions of human, mouse and rat BMP-4 share 98% aa sequence identity.

BMPs signal via different hetero-oligomeric complexes of type I and type II serine/threonine kinase receptors (for reviews, see references 9 - 10). BMP-4 receptors include the type I receptor, ALK-6/BMP RIB (11), and the type II receptor, BMP RII (12, 13). Signals from activated BMP receptors are directly transduced to the cell nucleus by Smad proteins that then become incorporated into transcriptional complexes (for a review, see reference 14). Smad1 and Smad4, for example, are components of the BMP-4-induced transcription complex that is essential for dorsoventral patterning in *Xenopus* embryos (15).

In addition to promoting bone formation, BMP-4 is involved in other aspects of development. BMP-4, produced by the dorsal aorta, can direct sympathetic neuron differentiation (16). It may also influence somite development by inhibiting the process of myogenesis (17). BMP-4 plays a central role in dorsal/ventral patterning (for a review, see reference 18). It specifies the development of ventral structures (*e.g.* skin from ectoderm and connective tissue/blood from mesoderm). The distribution of BMP-4 expression suggests a direct role in the specification of human hematopoietic cells from embryonic mesoderm *in vivo* (19). Dorsal structures (nervous systems and muscle) appear when BMP-4 signals are interrupted through the activities of binding proteins, such as Noggin.

Variable expression levels of BMP-4 have been linked to different pathological states. In mice, BMP-4 gene knockout by homologous recombination results in embryonic lethality (1). Cases of patients displaying an interstitial deletion of chromosome 14 suggest that the 14q22 region (*i.e.* BMP-4 gene location) is important for human eye and pituitary development (20). BMP-4 is overexpressed in patients suffering from fibrodysplasia ossificans progressiva (FOP, an extremely rare, inherited disorder that is associated with abnormal skeletal patterning) (21, 22).

The Quantikine BMP-4 Immunoassay is a 4.5 hour solid phase ELISA designed to measure BMP-4 levels in bone tissue extracts and cell culture supernates. It contains NS0-expressed, recombinant human BMP-4 and antibodies raised against the recombinant factor.

## PRINCIPLE OF THE ASSAY

This assay employs the quantitative sandwich enzyme immunoassay technique. A monoclonal antibody specific for BMP-4 has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and any BMP-4 present is bound by the immobilized antibody. After washing away any unbound substances, an enzyme-linked monoclonal antibody specific for BMP-4 is added to the wells. Following a wash to remove any unbound antibody-enzyme reagent, a substrate solution is added to the wells and color develops in proportion to the amount of BMP-4 bound. The color development is stopped and the intensity of the color is measured.

## LIMITATIONS OF THE PROCEDURE

- FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.
- The kit should not be used beyond the expiration date on the kit label.
- If samples generate values higher than the highest standard, further dilute the samples with Calibrator Diluent RD5-13 and repeat the assay.
- Any variation in standard diluent, operator, pipetting technique, washing technique, incubation time or temperature, and kit age can cause variation in binding.
- Cell lines tested did not have measureable levels of BMP-4 in this assay.
- This assay is designed to eliminate interference by soluble receptors, binding proteins, and other factor present in biological samples. Until all factors have been tested in this Immunoassay, the possibility of interference cannot be excluded.

## REAGENTS

**BMP-4 Microplate** (Part 890828) - 96 well polystyrene microplate (8 strips of 12 wells) coated with a monoclonal antibody against BMP-4.

**BMP-4 Conjugate** (Part 890829) - 21 mL of a monoclonal antibody against BMP-4 conjugated to horseradish peroxidase, with preservative.

**BMP-4 Standard** (Part 890830) - 20 ng of recombinant human BMP-4 in a buffered protein base with preservative, lyophilized.

**Assay Diluent RD1-61** (Part 895329) - 11 mL of a buffered protein base with preservative. May contain crystals. Bring to room temperature and mix gently to dissolve.

**Calibrator Diluent RD5-13** (Part 895309) - 21 mL of a buffered protein base with preservative.

**Wash Buffer Concentrate** (Part 895003) - 21 mL of a 25-fold concentrated solution of buffered surfactant with preservative.

**Color Reagent A** (Part 895000) - 12.5 mL of stabilized hydrogen peroxide.

**Color Reagent B** (Part 895001) - 12.5 mL of stabilized chromogen (tetramethylbenzidine).

**Stop Solution** (Part 895032) - 6 mL of 2 N sulfuric acid.

**Plate Covers** - 4 adhesive strips.

## STORAGE

<b>Unopened Kit</b>	Store at 2 - 8° C. Do not use past kit expiration date.	
<b>Opened/ Reconstituted Reagents</b>	Diluted Wash Buffer	May be stored for up to 1 month at 2 - 8° C.*
	Stop Solution	
	Assay Diluent RD1-61	
	Calibrator Diluent RD5-13	
	Conjugate	
	Unmixed Color Reagent A	
	Unmixed Color Reagent B	
	Standard	
	Microplate Wells	Return unused wells to the foil pouch containing the desiccant pack, reseal along entire edge of zip-seal. May be stored for up to 1 month at 2 - 8° C.*

\*Provided this is within the expiration date of the kit.

## OTHER SUPPLIES REQUIRED

- Microplate reader capable of measuring absorbance at 450 nm, with the correction wavelength set at 540 nm or 570 nm.
- Pipettes and pipette tips.
- Deionized or distilled water.
- Multi-channel pipette, squirt bottle, manifold dispenser, or automated microplate washer.
- 500 mL graduated cylinder.
- **12 mm x 75 mm polypropylene tubes.**

## PRECAUTION

The Stop Solution provided with this kit is an acid solution. Wear eye, hand, face, and clothing protection when using this material.

## SAMPLE COLLECTION AND STORAGE

**Warning: Polypropylene tubes must be used. Do not use glass.**

**Cell Culture Supernates** - Remove particulates by centrifugation and assay immediately or aliquot and store samples at  $\leq -20^{\circ}\text{C}$ . Avoid repeated freeze-thaw cycles.

**Bone** - Extract demineralized bone samples in 4 M Guanidine-HCl and protease inhibitors (23, 24). Dissolve the final sample in 2 M Guanidine-HCl.

**Note:** *Extractions can also be done in Urea (24, 25).*

Bone extract samples must be diluted in Calibrator Diluent RD5-13 prior to assay so that the concentration of Guanidine-HCl is  $\leq 0.5\text{ M}$ , and the concentration of Urea is  $\leq 1\text{ M}$ .

## REAGENT PREPARATION


**Bring all reagents to room temperature before use.**

**Wash Buffer** - If crystals have formed in the concentrate, warm to room temperature and mix gently until the crystals have completely dissolved. Dilute 20 mL of Wash Buffer Concentrate into deionized or distilled water to prepare 500 mL of Wash Buffer.

**Substrate Solution** - Color Reagents A and B should be mixed together in equal volumes within 15 minutes of use. Protect from light. 200  $\mu\text{L}$  of the resultant mixture is required per well.

**BMP-4 Standard** - Reconstitute the BMP-4 Standard with 1.0 mL of deionized water. This reconstitution produces a stock solution of 20,000 pg/mL. Allow the standard to sit for a minimum of 15 minutes with gentle agitation prior to making dilutions.

**Use polypropylene tubes.** Pipette 900  $\mu\text{L}$  of Calibrator Diluent RD5-13 into the 2000 pg/mL tube. Pipette 500  $\mu\text{L}$  of Calibrator Diluent RD5-13 into the remaining tubes. Use the stock solution to produce a dilution series (below). Mix each tube thoroughly before the next transfer. The 2000 pg/mL standard serves as the high standard. Calibrator Diluent RD5-13 serves as the zero standard (0 pg/mL).



## ASSAY PROCEDURE

**Bring all reagents and samples to room temperature before use. It is recommended that all samples and standards be assayed in duplicate.**

1. Prepare all reagents, working standards, and samples as directed in the previous sections.
2. Remove excess microplate strips from the plate frame, return them to the foil pouch containing the desiccant pack, and reseal tightly.
3. Add 100  $\mu\text{L}$  of Assay Diluent RD1-61 to each well. Assay Diluent RD1-61 may contain crystals. Warm to room temperature and mix gently to dissolve before use.
4. Add 50  $\mu\text{L}$  of Standard or sample\* per well. Cover with the adhesive strip provided. Incubate for 2 hours at room temperature.
5. Aspirate each well and wash, repeating the process three times for a total of four washes. Wash by filling each well with Wash Buffer (400  $\mu\text{L}$ ) using a squirt bottle, multi-channel pipette, manifold dispenser or autowasher. Complete removal of liquid at each step is essential to good performance. After the last wash, remove any remaining Wash Buffer by aspirating or by decanting. Invert the plate and blot it against clean paper towels.
6. Add 200  $\mu\text{L}$  of BMP-4 Conjugate to each well. Cover with a new adhesive strip. Incubate for 2 hours at room temperature. **Protect from light.**
7. Repeat the aspiration/wash as in step 5.
8. Add 200  $\mu\text{L}$  of Substrate Solution to each well. Incubate for 30 minutes at room temperature. **Protect from light.**
9. Add 50  $\mu\text{L}$  of Stop Solution to each well. If color change does not appear uniform, gently tap plate to ensure thorough mixing.
10. Determine the optical density of each well within 30 minutes, using a microplate reader set to 450 nm. If wavelength correction is available, set to 540 nm or 570 nm. If wavelength correction is not available, subtract readings at 540 nm or 570 nm from the readings at 450 nm. This subtraction will correct for optical imperfections in the plate. Readings made directly at 450 nm without correction may be higher and less accurate.

\*Bone samples require extraction and dilution. See Sample Collection section.

## ASSAY PROCEDURE SUMMARY

1. Prepare reagents, samples, and standards as instructed.



2. Add 100  $\mu\text{L}$  Assay Diluent RD1-61 to each well.  
Be sure crystals have dissolved before use.




3. Add 50  $\mu\text{L}$  Standard or sample\* to each well.  
Incubate 2 hours at RT.




4. Aspirate and wash 4 times.




5. Add 200  $\mu\text{L}$  Conjugate to each well.  
Incubate 2 hours at RT. **Protect from light.**



6. Aspirate and wash 4 times.



7. Add 200  $\mu\text{L}$  Substrate Solution to each well.  
Incubate 30 minutes at RT. **Protect from light.**



8. Add 50  $\mu\text{L}$  Stop Solution to each well.  
Read at 450 nm within 30 minutes.  
 $\lambda$  correction 540 or 570 nm

\*Bone samples require extraction and dilution.


## CALCULATION OF RESULTS

Average the duplicate readings for each standard, control, and sample and subtract the average zero standard optical density.

Create a standard curve by reducing the data using computer software capable of generating a four parameter logistic (4-PL) curve-fit. As an alternative, construct a standard curve by plotting the mean absorbance for each standard on the y-axis against the concentration on the x-axis and draw a best fit curve through the points on the graph. The data may be linearized by plotting the log of the BMP-4 concentrations versus the log of the O.D. and the best fit line can be determined by regression analysis. This procedure will produce an adequate but less precise fit of the data. If samples have been diluted, the concentration read from the standard curve must be multiplied by the dilution factor.

## TYPICAL DATA

This standard curve is provided for demonstration only. A standard curve should be generated for each set of samples assayed.



(pg/mL)	O.D.	Average	Corrected
0	0.008 0.008	0.008	-
31.2	0.072 0.065	0.068	0.060
62.5	0.140 0.128	0.134	0.126
125	0.263 0.243	0.253	0.245
250	0.454 0.449	0.452	0.444
500	0.851 0.793	0.822	0.814
1000	1.409 1.364	1.386	1.378
2000	2.366 2.435	2.400	2.392

## TECHNICAL HINTS

- Substrate Solution should remain colorless until added to the plate. Keep the Substrate Solution protected from light. Substrate Solution should change from colorless to gradations of blue. The color developed in the wells will turn from blue to yellow upon addition of the Stop Solution.
- Stop Solution should be added to the plate in the same order as the Substrate Solution.
- When mixing or reconstituting protein solutions, always avoid foaming.
- To avoid cross-contamination, change pipette tips between additions of each standard level, between sample additions, and between reagent additions. Also, use separate reservoirs for each reagent.
- To ensure accurate results, proper adhesion of plate sealers during incubation steps is necessary.

## PRECISION

### Intra-assay Precision (Precision within an assay)

Three samples of known concentration were tested twenty times on one plate to assess intra-assay precision.

### Inter-assay Precision (Precision between assays)

Three samples of known concentration were tested in forty separate assays to assess inter-assay precision.

Sample	Intra-assay Precision			Inter-assay Precision		
	1	2	3	1	2	3
n	20	20	20	40	40	40
Mean (pg/mL)	301	719	1351	289	734	1429
Standard deviation	15.9	22.5	60.1	16.9	42.7	75.4
CV (%)	5.3	3.1	4.4	5.8	5.8	5.3

## RECOVERY

The recovery of rhBMP-4 spiked to levels throughout the range of the assay in various matrices was evaluated.

Diluent	Average % Recovery	Range
2 M Guanidine-HCl* (n=2)	96	93 - 103%
2 M Urea* (n=2)	97	91 - 105%
Cell culture media (n=8)	100	91 - 109%

\*Samples were diluted prior to assay.

## SENSITIVITY

Thirty-six assays were evaluated and the minimum detectable dose (MDD) of BMP-4 ranged from 0.43 to 3.68 pg/mL. The mean MDD was 1.04 pg/mL.

The minimum detectable dose was determined by adding two standard deviations to the mean optical density value of twenty zero standard replicates and calculating the corresponding concentration.

## CALIBRATION

This immunoassay is calibrated against highly purified, NS0-expressed, recombinant human BMP-4 produced at R&D Systems.

## LINEARITY

The following samples were spiked with rhBMP-4 and then diluted with Calibrator Diluent RD5-13. The samples were then serially diluted and tested to assess linearity of the assay.

		Cell culture media (n=8)	2 M Guanidine-HCl* (n=2)	2 M Urea* (n=2)
1:2	Average % of Expected	104	102	101
	Range (%)	101-105	98-105	99-103
1:4	Average % of Expected	109	101	104
	Range (%)	105-111	100-103	101-107
1:8	Average % of Expected	110	101	98
	Range (%)	107-115	100-103	97-98
1:16	Average % of Expected	108	98	91
	Range (%)	103-111	94-101	90-92

## SPECIFICITY

This assay recognizes both natural and recombinant human BMP-4. The factors listed below were prepared at 50 ng/mL in Calibrator Diluent RD5-13 and assayed for cross-reactivity. Preparations of the following factors at 50 ng/mL in a mid-range recombinant human BMP-4 control were also assayed for interference. No significant cross-reactivity or interference was observed.

### Recombinant human:

Activin A  
Activin RI  
Activin RIIa  
Activin RIIb  
BMP-2  
BMP-5

BMP-6  
BMP-7  
BMPR-1A  
BMPR-1B  
IL-3 sR $\alpha$   
Inhibin A  
Inhibin B  
Lefty

Follistatin 288  
Follistatin 300  
Follistatin 315  
TGF- $\alpha$   
TGF- $\beta$ 1  
TGF- $\beta$ 2  
TGF- $\beta$ 3  
TGF- $\beta$ 5

TGF- $\beta$  sRI  
TGF- $\beta$  sRII  
**Recombinant mouse:**  
Follistatin  
SCF  
TNF- $\alpha$

### Recombinant rat:

Agrin  
**Other:**  
h $\alpha$ 2-macroglobulin  
hTGF- $\beta$ 1  
pTGF- $\beta$ 1

rmNoggin was found to interfere at concentrations > 1 ng/mL.

## REFERENCES

1. Reddi, A.H. (1994) *Curr. Opin. Genet. Dev.* **4**:737.
2. Reddi, A.H. (1998) *Nature Biotech.* **16**:247.
3. Leong, L.M. and P.M. Brickell (1996) *Int. J. Biochem. Cell Biol.* **28**:1293.
4. Wozney, J.M. (1992) *Mol. Reprod. Dev.* **32**:160.
5. Wozney, J.M. (1989) *Prog. Growth Factor Res.* **1**:267.
6. Mehler, M.F. *et al.* (1997) *Trends Neurosci.* **20**:309.
7. Sampath, T.K. *et al.* (1990) *J. Biol. Chem.* **265**:13198.
8. Wozney, J.M. *et al.* (1988) *Science* **242**:1528.
9. Miyazono, K. (1999) *Bone* **25**:91.
10. Yamashita, H. *et al.* (1996) *Bone* **19**:569.
11. ten Dijke, P. *et al.* (1994) *J. Biol. Chem.* **269**:16985.
12. Rosenzweig, B.L. *et al.* (1995) *Proc. Natl. Acad. Sci. USA* **92**:7632.
13. Nohno, T. *et al.* (1995) *J. Biol. Chem.* **270**:22522.
14. Kawabata, M. *et al.* (1998) *Cytokine Growth Factor Rev.* **9**:49.
15. Henningfeld, K.A. *et al.* (2000) *J. Biol. Chem.* **275**:21827.
16. Reissman, E. *et al.* (1996) *Development* **122**:2079.
17. Tajbakhsh, S. and G. Cossu (1997) *Curr. Opin. Genet. Dev.* **7**:634.
18. Graff, J.M. (1997) *Cell* **89**:171.
19. Marshall, C.J. *et al.* (2000) *Blood* **96**:1591.
20. Lemyre, E. *et al.* (1998) *Am. J. Med. Genet.* **77**:162.
21. Xu, M. and E.M. Shore (1998) *Clin. Orthop.* **346**:53.
22. Shore, E.M. *et al.* (1998) *Calcif. Tissue Int.* **63**:221.
23. Takaoka, K. *et al.*, (1980) *Clin. Orthop. Relat. Res.* **148**:274.
24. Sampath, T.K. and A.H. Reddi, (1981) *Proc. Natl. Acad. Sci. USA* **78**(12):7599.
25. Urist, M.R. *et al.*, (1982) *Clin. Orthop. Relat. Res.* **162**:219.