

Quantikine[®] ELISA

Human IFN- γ Immunoassay

Catalog Number DIF50
SIF50
PDIF50

For the quantitative determination of human interferon gamma (IFN- γ) concentrations in cell culture supernates, serum, and plasma.

This package insert must be read in its entirety before using this product.
For research use only. Not for use in diagnostic procedures.

TABLE OF CONTENTS

| SECTION | PAGE |
|---|------|
| INTRODUCTION | 1 |
| PRINCIPLE OF THE ASSAY | 2 |
| LIMITATIONS OF THE PROCEDURE | 2 |
| TECHNICAL HINTS | 2 |
| MATERIALS PROVIDED & STORAGE CONDITIONS | 3 |
| OTHER SUPPLIES REQUIRED | 4 |
| PRECAUTION | 4 |
| SAMPLE COLLECTION & STORAGE | 4 |
| REAGENT PREPARATION | 5 |
| ASSAY PROCEDURE | 6 |
| CALCULATION OF RESULTS | 7 |
| TYPICAL DATA | 7 |
| PRECISION | 8 |
| RECOVERY | 8 |
| SENSITIVITY | 8 |
| LINEARITY | 9 |
| CALIBRATION | 9 |
| SAMPLE VALUES | 9 |
| SPECIFICITY | 10 |
| REFERENCES | 11 |
| PLATE LAYOUT | 12 |

MANUFACTURED AND DISTRIBUTED BY:

USA & Canada | R&D Systems, Inc.

614 McKinley Place NE, Minneapolis, MN 55413, USA
TEL: (800) 343-7475 (612) 379-2956 FAX: (612) 656-4400
E-MAIL: info@RnDSystems.com

DISTRIBUTED BY:

UK & Europe | R&D Systems Europe, Ltd.

19 Barton Lane, Abingdon Science Park, Abingdon OX14 3NB, UK
TEL: +44 (0)1235 529449 FAX: +44 (0)1235 533420
E-MAIL: info@RnDSystems.co.uk

China | R&D Systems China Co., Ltd.

24A1 Hua Min Empire Plaza, 726 West Yan An Road, Shanghai PRC 200050
TEL: +86 (21) 52380373 FAX: +86 (21) 52371001
E-MAIL: info@RnDSystemsChina.com.cn

INTRODUCTION

Interferon-gamma (IFN- γ , also known as type II interferon) is an important immunoregulatory cytokine that was originally identified through its anti-viral activity (1, 2). It plays key roles in host defense by exerting anti-viral, anti-proliferative, and immunoregulatory activities (3, 4). On many cell types, IFN- γ induces the production of cytokines and upregulates the expression of various membrane proteins including class I and II MHC antigens, Fc receptors, leukocyte adhesion molecules, and B7 family antigens. IFN- γ is a potent activator of macrophage effector functions. It directs the synthesis, class switching, and secretion of immunoglobulins by B cells. IFN- γ also influences T-helper cell phenotype development by inhibiting Th2 differentiation and stimulating Th1 development (3, 4). IFN- γ plays a central role in the progression of inflammatory diseases such as autoimmunity and atherosclerosis (5, 6).

Biologically active IFN- γ consists of a noncovalently linked homodimer of 20-25 kDa variably glycosylated subunits (7). Mature human IFN- γ shares 90% amino acid (aa) sequence identity with rhesus IFN- γ , 59-64% with bovine, canine, equine, feline, and porcine IFN- γ , and 37-43% aa identity with cotton rat, mouse, and rat IFN- γ . IFN- γ dimers bind to transmembrane IFN- γ RI (alpha subunits) which then interact with transmembrane IFN- γ RII (beta subunits) to form a functional receptor complex of two α and two β subunits (8, 9). Inclusion of IFN- γ RII in the receptor complex increases the ligand binding affinity as well as the efficiency of signal transduction (9, 10). Whereas the α -chain is expressed constitutively on many cell types, the cellular regulation of the β -chain correlates with an IFN- γ responsive state and is tightly regulated (8).

IFN- γ is produced by a number of cell types under inflammatory conditions, including dendritic epidermal/ $\gamma\delta$ T cells (11), keratinocytes (12), peripheral blood $\gamma\delta$ T cells (13), mast cells (14), neurons (15), CD8⁺ T cells (16), macrophages (17), B cells (18), neutrophils (19), NK cells (20), CD4⁺ T cells (21), and testicular spermatids (22).

The Quantikine Human IFN- γ Immunoassay is a 4.5 hour solid phase ELISA designed to measure IFN- γ levels in cell culture supernates, serum, and plasma. It contains *E. coli*-expressed recombinant human IFN- γ and antibodies raised against the recombinant factor. Results obtained for naturally occurring human IFN- γ samples showed linear curves that were parallel to the standard curves obtained using the Quantikine kit standards. These results indicate that this kit can be used to determine relative mass values for natural human IFN- γ .

The presence of a naturally occurring form of a receptor for IFN- γ has been reported in normal human urine (23). Until the soluble IFN- γ receptor has been tested in the Quantikine Human IFN- γ Immunoassay, the possibility of interference cannot be excluded.

PRINCIPLE OF THE ASSAY

This assay employs the quantitative sandwich enzyme immunoassay technique. Polyclonal antibody specific for IFN- γ has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and any IFN- γ present is bound by the immobilized antibody. After washing away any unbound substances, an enzyme-linked polyclonal antibody specific for IFN- γ is added to the wells. Following a wash to remove any unbound antibody-enzyme reagent, a substrate solution is added to the wells and color develops in proportion to the amount of IFN- γ bound in the initial step. The color development is stopped and the intensity of the color is measured.

LIMITATIONS OF THE PROCEDURE

- FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.
- The kit should not be used beyond the expiration date on the kit label.
- Do not mix or substitute reagents with those from other lots or sources.
- If samples generate values higher than the highest standard, dilute the samples with Calibrator Diluent and repeat the assay. If cell culture supernate samples require a large dilution, perform an intermediate dilution with culture media and the final dilution in Calibrator Diluent.
- Any variation in standard diluent, operator, pipetting technique, washing technique, incubation time or temperature, and kit age can cause variation in binding.
- Variations in sample collection, processing, and storage may cause sample value differences.
- This assay is designed to eliminate interference by soluble receptors, binding proteins, and other factors present in biological samples. Until all factors have been tested in the Quantikine Immunoassay, the possibility of interference cannot be excluded.

TECHNICAL HINTS

- When mixing or reconstituting protein solutions, always avoid foaming.
- To avoid cross-contamination, change pipette tips between additions of each standard level, between sample additions, and between reagent additions. Also, use separate reservoirs for each reagent.
- To ensure accurate results, proper adhesion of plate sealers during incubation steps is necessary.
- When using an automated plate washer, adding a 30 second soak period following the addition of Wash Buffer, and/or rotating the plate 180 degrees between wash steps may improve assay precision.
- Substrate Solution should remain colorless until added to the plate. Keep Substrate Solution protected from light. Substrate Solution should change from colorless to gradations of blue.
- Stop Solution should be added to the plate in the same order as the Substrate Solution. The color developed in the wells will turn from blue to yellow upon addition of the Stop Solution. Wells that are green in color indicate that the Stop Solution has not mixed thoroughly with the Substrate Solution.

MATERIALS PROVIDED & STORAGE CONDITIONS

Store the unopened kit at 2-8 °C. Do not use past kit expiration date.

| PART | PART # | CATALOG # DIF50 | CATALOG # SIF50 | DESCRIPTION | STORAGE OF OPENED/ RECONSTITUTED MATERIAL |
|------------------------------|--------|--------------------|--------------------|--|--|
| IFN- γ Microplate | 890582 | 1 plate | 6 plates | 96 well polystyrene microplate (12 strips of 8 wells) coated with polyclonal antibody against IFN- γ . | Return unused wells to the foil pouch containing the desiccant pack. Reseal along entire edge of zip-seal. May be stored for up to 1 month at 2-8 °C.* |
| IFN- γ Standard | 890210 | 1 vial | 6 vials | Recombinant human IFN- γ in a buffered protein base with preservatives; lyophilized. | Aliquot and store at ≤ -20 °C for up to 1 month in a manual defrost freezer.* Avoid repeated freeze-thaw cycles. |
| IFN- γ Conjugate | 890583 | 1 vial | 6 vials | 21 mL/vial of polyclonal antibody against IFN- γ conjugated to horseradish peroxidase with preservatives. | May be stored for up to 1 month at 2-8 °C.* |
| Assay Diluent RD1-51 | 895342 | 1 vial | 6 vials | 11 mL/vial of a buffered protein base with preservatives and blue dye. | |
| Calibrator Diluent RD6-21 | 895261 | 1 vial | 6 vials | 21 mL/vial of a buffered protein base with preservatives. | |
| Wash Buffer Concentrate | 895003 | 1 vial | 6 vials | 21 mL/vial of a 25-fold concentrated solution of buffered surfactant with preservatives. | |
| Color Reagent A | 895000 | 1 vial | 6 vials | 12 mL/vial of stabilized hydrogen peroxide. | |
| Color Reagent B | 895001 | 1 vial | 6 vials | 12 mL/vial of stabilized chromogen (tetramethylbenzidine). | |
| Stop Solution | 895032 | 1 vial | 6 vials | 6 mL/vial of 2 N sulfuric acid. | |
| Plate Sealers | N/A | 4 strips | 24 strips | Adhesive strips. | |

* Provided this is within the expiration date of the kit.

DIF50 contains sufficient materials to run an ELISA on one 96 well plate.

SIF50 (SixPak) contains sufficient materials to run ELISAs on six 96 well plates.

This kit is also available in a PharmPak (R&D Systems, Catalog # PDIF50). PharmPaks contain sufficient materials to run ELISAs on 50 microplates. Specific vial counts of each component may vary. Please refer to the literature accompanying your order for specific vial counts.

OTHER SUPPLIES REQUIRED

- Microplate reader capable of measuring absorbance at 450 nm, with the correction wavelength set at 540 nm or 570 nm.
- Pipettes and pipette tips.
- 500 mL graduated cylinder.
- Deionized or distilled water.
- Squirt bottle, manifold dispenser, or automated microplate washer.
- **Polypropylene** test tubes for dilution of standards.
- Human IFN- γ controls (optional; available from R&D Systems).

PRECAUTION

The Stop Solution provided with this kit is an acid solution. Wear protective gloves, clothing, eye, and face protection. Wash hands thoroughly after handling.

SAMPLE COLLECTION & STORAGE

The sample collection and storage conditions listed below are intended as general guidelines. Sample stability has not been evaluated.

Cell Culture Supernates - Remove particulates by centrifugation and assay immediately or aliquot and store samples at ≤ -20 °C. Avoid repeated freeze-thaw cycles.

Serum - Use a serum separator tube (SST) and allow samples to clot for 30 minutes before centrifugation for 15 minutes at 1000 x g. Remove serum and assay immediately or aliquot and store samples at ≤ -20 °C. Avoid repeated freeze-thaw cycles.

Plasma - Collect plasma using EDTA as an anticoagulant. Centrifuge for 15 minutes at 1000 x g within 30 minutes of collection. Assay immediately or aliquot and store samples at ≤ -20 °C. Avoid repeated freeze-thaw cycles. **Heparin and citrate plasma samples cannot be used in this assay.**

Note: *Hemolyzed samples are not suitable for the measurement of human IFN- γ with this assay. Samples with high levels of rheumatoid factor may interfere with the measurement of human IFN- γ in this assay.*

REAGENT PREPARATION


Bring all reagents to room temperature before use.

Wash Buffer - If crystals have formed in the concentrate, warm to room temperature and mix gently until the crystals have completely dissolved. Dilute 20 mL of Wash Buffer Concentrate into deionized or distilled water to prepare 500 mL of Wash Buffer.

Substrate Solution - Color Reagents A and B should be mixed together in equal volumes within 15 minutes of use. Protect from light. 200 μL of the resultant mixture is required per well.

IFN- γ Standard - **Refer to the vial label for reconstitution volume.** Reconstitute the IFN- γ Standard with Calibrator Diluent RD6-21. This reconstitution produces a stock solution of 1000 pg/mL. Allow the standard to sit for a minimum of 15 minutes with gentle agitation prior to making dilutions.

Use polypropylene tubes. Pipette 500 μL of Calibrator Diluent RD6-21 into each tube. Use the stock solution to produce a dilution series (below). Mix each tube thoroughly before the next transfer. The undiluted standard serves as the high standard (1000 pg/mL). Calibrator Diluent RD6-21 serves as the zero standard (0 pg/mL).



ASSAY PROCEDURE

Bring all reagents and samples to room temperature before use. It is recommended that all standards, samples, and controls be assayed in duplicate.

1. Prepare all reagents and working standards as directed in the previous sections.
2. Remove excess microplate strips from the plate frame, return them to the foil pouch containing the desiccant pack, and reseal.
3. Add 100 μL of Assay Diluent RD1-51 to each well.
4. Add 100 μL of Standard, sample, or control per well. **Ensure reagent addition is uninterrupted and completed within 15 minutes.** Cover with the adhesive strip provided. Incubate for 2 hours at room temperature. A plate layout is provided to record standards and samples assayed.
5. Aspirate each well and wash, repeating the process three times for a total of four washes. Wash by filling each well with Wash Buffer (400 μL) using a squirt bottle, manifold dispenser, or autowasher. Complete removal of liquid at each step is essential to good performance. After the last wash, remove any remaining Wash Buffer by aspirating or decanting. Invert the plate and blot it against clean paper towels.
6. Add 200 μL of IFN- γ Conjugate to each well. Cover with a new adhesive strip. Incubate for 2 hours at room temperature.
7. Repeat the aspiration/wash as in step 5.
8. Add 200 μL of Substrate Solution to each well. Incubate for 30 minutes at room temperature. **Protect from light.**
9. Add 50 μL of Stop Solution to each well. The color in the well should change from blue to yellow. If the color in the well is green or if the color change does not appear uniform, gently tap the plate to ensure thorough mixing.
10. Determine the optical density of each well within 30 minutes, using a microplate reader set to 450 nm. If wavelength correction is available, set to 540 nm or 570 nm. If wavelength correction is not available, subtract readings at 540 nm or 570 nm from the readings at 450 nm. This subtraction will correct for optical imperfections in the plate. Readings made directly at 450 nm without correction may be higher and less accurate.

CALCULATION OF RESULTS


Average the duplicate readings for each standard, control, and sample and subtract the average zero standard optical density.

Create a standard curve by reducing the data using computer software capable of generating a four parameter logistic (4-PL) curve-fit. As an alternative, construct a standard curve by plotting the mean absorbance for each standard on the y-axis against the concentration on the x-axis and draw a best fit curve through the points on the graph. The data may be linearized by plotting the log of the IFN- γ concentrations versus the log of the O.D. and the best fit line can be determined by regression analysis. This procedure will produce an adequate but less precise fit of the data.

If samples have been diluted, the concentration read from the standard curve must be multiplied by the dilution factor.

TYPICAL DATA

This standard curve is provided for demonstration only. A standard curve should be generated for each set of samples assayed.



| (pg/mL) | O.D. | Average | Corrected |
|---------|----------------|---------|-----------|
| 0 | 0.094 0.093 | 0.094 | — |
| 15.6 | 0.125 0.129 | 0.127 | 0.033 |
| 31.2 | 0.163 0.162 | 0.162 | 0.068 |
| 62.5 | 0.231 0.235 | 0.233 | 0.139 |
| 125 | 0.373 0.376 | 0.374 | 0.280 |
| 250 | 0.644 0.670 | 0.657 | 0.563 |
| 500 | 1.190 1.187 | 1.188 | 1.094 |
| 1000 | 2.196 2.237 | 2.216 | 2.122 |

PRECISION

Intra-assay Precision (Precision within an assay)

Three samples of known concentration were tested twenty times on one plate to assess intra-assay precision.

Inter-assay Precision (Precision between assays)

Three samples of known concentration were tested in forty separate assays to assess inter-assay precision.

| Sample | Intra-Assay Precision | | | Inter-Assay Precision | | |
|--------------------|-----------------------|-----|------|-----------------------|------|------|
| | 1 | 2 | 3 | 1 | 2 | 3 |
| n | 20 | 20 | 20 | 40 | 40 | 40 |
| Mean (pg/mL) | 79.2 | 214 | 466 | 78.7 | 206 | 457 |
| Standard deviation | 3.7 | 5.5 | 13.0 | 6.1 | 13.1 | 17.0 |
| CV (%) | 4.7 | 2.6 | 2.8 | 7.8 | 6.4 | 3.7 |

RECOVERY

The recovery of IFN- γ spiked to levels throughout the range of the assay in various matrices was evaluated.

| Sample Type | Average % Recovery | Range |
|--------------------------|--------------------|---------|
| Cell culture media (n=4) | 102 | 93-109% |
| Serum (n=5) | 102 | 91-118% |
| EDTA plasma (n=5) | 98 | 88-111% |

SENSITIVITY

The minimum detectable dose (MDD) of IFN- γ is typically less than 8.0 pg/mL.

The MDD was determined by adding two standard deviations to the mean optical density value of twenty zero standard replicates and calculating the corresponding concentration.

LINEARITY

To assess the linearity of the assay, samples were spiked with high concentrations of IFN- γ in various matrices and diluted with Calibrator Diluent RD6-21 to produce samples with values within the dynamic range of the assay.

| | | Cell culture media (n=5) | Serum (n=5) | EDTA plasma (n=5) |
|------|-----------------------|--------------------------|-------------|-------------------|
| 1:2 | Average % of Expected | 98 | 100 | 104 |
| | Range (%) | 95-103 | 98-101 | 101-109 |
| 1:4 | Average % of Expected | 100 | 99 | 107 |
| | Range (%) | 93-106 | 97-103 | 101-115 |
| 1:8 | Average % of Expected | 98 | 96 | 101 |
| | Range (%) | 88-109 | 92-99 | 94-107 |
| 1:16 | Average % of Expected | 93 | 89 | 94 |
| | Range (%) | 75-109 | 84-98 | 86-99 |

CALIBRATION

This immunoassay is calibrated against a highly purified *E. coli*-expressed recombinant human IFN- γ produced at R&D Systems.

The NIBSC/WHO 1st British Standard for human leukocyte IFN- γ (82/587) was evaluated in this kit in July 2008. To convert sample values obtained with the Quantikine Human IFN- γ kit to approximate NIBSC units, use the equation below:

NIBSC/WHO (82/587) approximate value (U/mL) = 0.017 x Quantikine Human IFN- γ value (pg/mL)

SAMPLE VALUES

Serum/Plasma - Thirty samples from apparently healthy volunteers were evaluated for the presence of IFN- γ in this assay. No medical histories were available for the donors used in this study. All samples measured less than the lowest IFN- γ standard, 15.6 pg/mL.

Cell Culture Supernates - Human peripheral blood mononuclear cells (1×10^6 cells/mL) were cultured in RPMI supplemented with 10% fetal calf serum, 50 μ M β -mercaptoethanol, 2 mM L-glutamine, 100 U/mL penicillin, and 100 μ g/mL streptomycin sulfate. Cells were cultured unstimulated or stimulated with 10 μ g/mL of PHA. Aliquots of the cell culture supernates were removed on days 1 and 5 and assayed for levels of natural IFN- γ .

| Condition | Day 1 (pg/mL) | Day 5 (pg/mL) |
|--------------|---------------|---------------|
| Unstimulated | 453 | 2651 |
| Stimulated | 8360 | 7851 |

SPECIFICITY

This assay recognizes natural and recombinant human IFN- γ .

The factors listed below were prepared at 50 ng/mL in Calibrator Diluent RD6-21 and assayed for cross-reactivity. Preparations of the following factors at 50 ng/mL in a mid-range rhIFN- γ control were assayed for interference. No significant cross-reactivity or interference was observed.

Recombinant human:

IFN- γ R1

Recombinant mouse:

IFN- γ

IFN- γ R1

Other recombinants:

bovine IFN- γ

canine IFN- γ

cotton rat IFN- γ

equine IFN- γ

feline IFN- γ

porcine IFN- γ

rat IFN- γ

rhesus macaque IFN- γ

Recombinant porcine IFN- γ was found to cross-react approximately 0.08%.


Recombinant rhesus macaque IFN- γ was found to cross-react approximately 0.12%.

REFERENCES

1. Billiau, A. and P. Matthys (2009) Cytokine Growth Factor Rev. **20**:97.
2. Wheelock, E.F. (1965) Science **146**:310.
3. Schoenborn, J.R and C.B. Wilson (2007) Adv. Immunol. **96**:41.
4. Pestka, S. *et al.* (2004) Immunol. Rev. **202**:8.
5. Kelchtermans, H. *et al.* (2008) Trends Immunol. **29**:479.
6. McLaren, J.E. and D.P. Ramji (2009) Cytokine Growth Factor Rev. **20**:125.
7. Gray, P.W. and D.V. Goeddel (1982) Nature **298**:859.
8. Bach, E.A. *et al.* (1997) Annu. Rev. Immunol. **15**:563.
9. Marsters, S.A. *et al.* (1995) Proc. Natl. Acad. Sci. **92**:5401.
10. Krause, C.D. *et al.* (2000) J. Biol. Chem. **275**:22995.
11. Sugaya, M. *et al.* (1999) J. Invest. Dermatol. **113**:350.
12. Howie, S.E.M. *et al.* (1996) J. Invest. Dermatol. **106**:1218.
13. Battistini, L. *et al.* (1997) J. Immunol. **159**:3723.
14. Gupta, A.A. *et al.* (1996) J. Immunol. **157**:2123.
15. Neumann, H. *et al.* (1997) J. Exp. Med. **186**:2023.
16. Hoiden, I. and G. Moller (1996) Scand. J. Immunol. **44**:501.
17. Puddu, P. *et al.* (1997) J. Immunol. **159**:3490.
18. Yoshimoto, T. *et al.* (1997) Proc. Natl. Acad. Sci. USA **94**:3948.
19. Yeaman, G.R. *et al.* (1998) J. Immunol. **160**:5145.
20. Asea, A. *et al.* (1996) Clin. Exp. Immunol. **105**:376.
21. Briscoe, D.M. *et al.* (1997) J. Immunol. **159**:3247.
22. Dejuco, N. *et al.* (1995) Endocrinology **136**:4925.
23. Novick, D. *et al.* (1989) J. Exp. Med. **170**:1409.

PLATE LAYOUT

Use this plate layout to record standards and samples assayed.



NOTES

NOTES