

DESCRIPTION

Source *E. coli*-derived
Phe25-Leu211, with an N-terminal Met
Accession # P08505

N-terminal Sequence Analysis Phe25

Predicted Molecular Mass 21.8 kDa

SPECIFICATIONS

Activity Measured in a cell proliferation assay using T1165.85.2.1 mouse plasmacytoma cells. Nordan, R.P. *et al.* (1987) *J. Immunol.* **139**:813. The ED₅₀ for this effect is typically 0.02-0.06 ng/mL.

Endotoxin Level <1.0 EU per 1 µg of the protein by the LAL method.

Purity >97%, by SDS-PAGE under reducing conditions and visualized by silver stain.

Formulation Lyophilized from a 0.2 µm filtered solution in Sodium Acetate and EDTA with Trehalose and with BSA as a carrier protein. See Certificate of Analysis for details.

PREPARATION AND STORAGE

Reconstitution Reconstitute at 100 µg/mL in sterile PBS containing at least 0.1% human or bovine serum albumin.

Shipping The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature recommended below.

Stability & Storage **Use a manual defrost freezer and avoid repeated freeze-thaw cycles.**

- 12 months from date of receipt, -20 to -70 °C as supplied.
- 1 month, 2 to 8 °C under sterile conditions after reconstitution.
- 3 months, -20 to -70 °C under sterile conditions after reconstitution.

BACKGROUND

Interleukin 6 (IL-6) is a pleiotropic α-helical cytokine that plays important roles in acute phase reactions, inflammation, hematopoiesis, bone metabolism, and cancer progression. IL-6 activity is central to the transition from acute inflammation to either acquired immunity or chronic inflammatory disease. It is secreted by multiple cell types as a 22 kDa – 28 kDa phosphorylated and variably glycosylated molecule (1 – 4). Mature mouse IL-6 is 187 amino acids (aa) in length and shares 42% and 85% aa sequence identity with human and rat IL-6, respectively (5). Alternate splicing generates several isoforms with internal deletions (6). Mouse IL-6 is equally active on human and rat cells (7). IL-6 induces signaling through a cell surface heterodimeric receptor complex composed of a ligand binding subunit (IL-6 R) and a signal transducing subunit (gp130). IL-6 binds to IL-6 R, triggering IL-6 R association with gp130 and gp130 dimerization (8). gp130 is also a component of the receptors for CLC, CNTF, CT-1, IL-11, IL-27, LIF, and OSM (9). Soluble forms of IL-6 R are generated by both alternate splicing and proteolytic cleavage (9). In a mechanism known as trans-signaling, complexes of soluble IL-6 and IL-6 R elicit responses from gp130-expressing cells that lack cell surface IL-6 R (3). Trans-signaling enables a wider range of cell types to respond to IL-6, as the expression of gp130 is ubiquitous while that of IL-6 R is predominantly restricted to hepatocytes, leukocytes, and lymphocytes (3). Soluble splice forms of gp130 block trans-signaling from IL-6/IL-6 R but not from other cytokines that utilize gp130 as a coreceptor (4, 10).

References:

1. Van Snick, J. (1990) *Annu. Rev. Immunol.* **8**:253.
2. Hodge, D.R. *et al.* (2005) *Eur. J. Cancer* **41**:2502.
3. Jones, S.A. (2005) *J. Immunol.* **175**:3468.
4. Rose-John, S. *et al.* (2006) *J. Leukoc. Biol.* **80**:227.
5. Van Snick, J. *et al.* (1988) *Eur. J. Immunol.* **18**:193.
6. Yatsenko, O.P. *et al.* (2004) *Cytokine* **28**:190.
7. Chiu, C.P. *et al.* (1988) *Proc. Natl. Acad. Sci.* **85**:7099.
8. Murakami, M. *et al.* (1993) *Science* **260**:1808.
9. Muller-Newen, G. (2003) *Sci. STKE* **2003**:PE40.
10. Mitsuyama, K. *et al.* (2006) *Clin. Exp. Immunol.* **143**:125.