

DESCRIPTION

Source	<i>E. coli</i> -derived Leu80-Leu235, with an N-terminal Met Accession # P06804
N-terminal Sequence Analysis	Met & Ser84
Structure / Form	A small amount of recombinant protein lacking the N-terminal methionine and four additional amino acid residues may be also present.
Predicted Molecular Mass	17 kDa

SPECIFICATIONS

Activity	Measured in a cytotoxicity assay using L-929 mouse fibroblast cells in the presence of the metabolic inhibitor actinomycin D. Matthews, N. and M.L. Neale (1987) in <i>Lymphokines and Interferons, A Practical Approach</i> . Clemens, M.J. <i>et al.</i> (eds): IRL Press. 221. The ED ₅₀ for this effect is typically 0.008-0.05 ng/mL.
Endotoxin Level	<1.0 EU per 1 μ g of the protein by the LAL method.
Purity	>97%, by SDS-PAGE under reducing conditions and visualized by silver stain.
Formulation	Lyophilized from a 0.2 μ m filtered solution in PBS with BSA as a carrier protein. See Certificate of Analysis for details.

PREPARATION AND STORAGE

Reconstitution	Reconstitute at 100 μ g/mL in sterile PBS containing at least 0.1% human or bovine serum albumin.
Shipping	The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature recommended below.
Stability & Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. <ul style="list-style-type: none"> ● 12 months from date of receipt, -20 to -70 °C as supplied. ● 1 month, 2 to 8 °C under sterile conditions after reconstitution. ● 3 months, -20 to -70 °C under sterile conditions after reconstitution.

BACKGROUND

Tumor necrosis factor alpha (TNF- α , also known as cachectin and TNFSF1A), is the prototypic ligand of the TNF superfamily. It is a pleiotropic molecule that plays a central role in inflammation, apoptosis, and immune system development. TNF- α is produced by a wide variety of immune and epithelial cell types (1, 2). Mouse TNF- α consists of a 35 amino acid (aa) cytoplasmic domain, a 21 aa transmembrane segment, and a 179 aa extracellular domain (ECD) (3). Within the ECD, mouse TNF- α shares 94% aa sequence identity with rat and 70% - 77% with bovine, canine, cotton rat, equine, feline, human, porcine, rat, and rhesus TNF- α . The 26 kDa type 2 transmembrane protein is assembled intracellularly to form a noncovalently linked homotrimer (4). Ligation of this complex induces reverse signaling that promotes lymphocyte costimulation but diminishes monocyte responsiveness (5). Cleavage of membrane bound TNF- α by TACE/ADAM17 releases a 55 kDa soluble trimeric form of TNF- α (6, 7). TNF- α trimers bind the ubiquitous TNF RI and the hematopoietic cell-restricted TNF RII, both of which are also expressed as homotrimers (1, 8). TNF- α regulates lymphoid tissue development through control of apoptosis (2). It also promotes inflammatory responses by inducing the activation of vascular endothelial cells and macrophages (2). TNF- α is a key cytokine in the development of several inflammatory disorders (9). It contributes to the development of type 2 diabetes through its effects on insulin resistance and fatty acid metabolism (10, 11).

References:

1. Idriss, H.T. and J.H. Naismith (2000) *Microsc. Res. Tech.* **50**:184.
2. Hehlhans, T. and K. Pfeffer (2005) *Immunology* **115**:1.
3. Fransen, L. *et al.* (1985) *Nucl. Acids Res.* **13**:4417.
4. Tang, P. *et al.* (1996) *Biochemistry* **35**:8216.
5. Eissner G. *et al.* (2004) *Cytokine Growth Factor Rev.* **15**:353.
6. Black, R.A. *et al.* (1997) *Nature* **385**:729.
7. Moss, M.L. *et al.* (1997) *Nature* **385**:733.
8. Loetscher, H. *et al.* (1991) *J. Biol. Chem.* **266**:18324.
9. Clark, I.A. (2007) *Cytokine Growth Factor Rev.* **18**:335.
10. Romanatto, T. *et al.* (2007) *Peptides* **28**:1050.
11. Hector, J. *et al.* (2007) *Horm. Metab. Res.* **39**:250.