

Biotinylated Anti-human SHBG Antibody

ORDERING INFORMATION

Catalog Number: BAF2656

Lot Number: VQO01

Size: 50 µg

Formulation: 0.2 µm filtered solution in PBS with BSA

Storage: -20° C

Reconstitution: sterile 0.1% BSA in TBS

Specificity: human SHBG

Immunogen: NS0-derived rhSHBG

Ig Type: goat IgG

Applications: Western blot
Immunohistochemistry

Preparation

Produced in goats immunized with purified, NS0-derived, recombinant human Sex Hormone Binding Globulin (rhSHBG). Human SHBG specific IgG was purified by human SHBG affinity chr and then biotinylated. Human SHBG (also known as ABP or androgen-binding protein) is a variably glycosylated, secreted non-disulfide linked homodimer that belongs to the SHBG family. Members of this small family have tandem repeats of 170 amino acid (aa) long laminin α chain G-like domains. Each SHBG monomer is a 47 - 53 kDa, 373 aa glycoprotein that contains one steroid-binding site in its N-terminal G-like domain. Male and female sex hormones are bound with equal affinity. SHBG is synthesized by the liver and circulates in blood; ABP (SHBG in the testis) is synthesized by Sertoli cells and circulates in the male reproductive system. Human SHBG shares approximately 67% aa sequence identity with mouse and rat SHBG.

Formulation

Lyophilized from a 0.2 µm filtered solution in phosphate-buffered saline (PBS) containing 50 µg of bovine serum albumin (BSA) per 1 µg of antibody.

Reconstitution

Reconstitute with sterile Tris-buffered saline pH 7.3 (20 mM Trizma base, 150 mM NaCl) containing 0.1% BSA. If 1 mL of buffer is used, the antibody concentration will be 50 µg/mL.

Storage

Lyophilized samples are stable for twelve months from date of receipt when stored at -20° C to -70° C. Upon reconstitution, the antibody can be stored at 2° - 8° C for 1 month without detectable loss of activity. Reconstituted antibody can also be aliquotted and stored frozen at -20° C to -70° C **in a manual defrost freezer** for six months without detectable loss of activity. **Avoid repeated freeze-thaw cycles.**

Specificity

This antibody has been selected for use as a detection antibody in the applications listed below.

Applications

Western Blot - This antibody can be used at 0.1 - 0.2 µg/mL with the appropriate secondary reagents to detect human SHBG. The detection limit for rhSHBG is approximately 5 ng/lane under non-reducing and reducing conditions.

Immunohistochemistry - This antibody will detect SHBG in cells and tissues. The working dilution is 2 - 15 µg/mL. For chromogenic detection of labeling, use R&D Systems' Cell and Tissue Staining Kits (CTS Series).

Optimal dilutions should be determined by each laboratory for each application.