

DESCRIPTION

Species Reactivity	Human
Specificity	Detects recombinant human Fetuin A in direct ELISAs and Western blots. Does not cross-react with 16 other related proteins tested.
Source	Monoclonal Mouse IgG ₁ Clone # 162922
Purification	Protein A or G purified from hybridoma culture supernatant
Immunogen	Mouse myeloma cell line NS0-derived recombinant human Fetuin A Ala19 - Val367 Accession # P02765
Formulation	Lyophilized from a 0.2 µm filtered solution in PBS with Trehalose. See Certificate of Analysis for details.

APPLICATIONS

Please Note: Optimal dilutions should be determined by each laboratory for each application. [General Protocols](#) are available in the *Technical Information* section on our website.

	Recommended Concentration	Sample
Western Blot	1 µg/mL	Recombinant Human Fetuin A/AHSG (Catalog # 1184-PI)

PREPARATION AND STORAGE

Reconstitution	Reconstitute at 0.5 mg/mL in sterile PBS.
Shipping	The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature recommended below.
Stability & Storage	<p>Use a manual defrost freezer and avoid repeated freeze-thaw cycles.</p> <ul style="list-style-type: none"> ● 12 months from date of receipt, -20 to -70 °C as supplied. ● 1 month from date of receipt, 2 to 8 °C, reconstituted. ● 6 months from date of receipt, -20 to -70 °C, reconstituted.

BACKGROUND

Fetuin A, also named Alpha 2-Heremans Schmid Glycoprotein (AHSG), is an abundant plasma protein that is a member of the cystatin superfamily.